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Dear Customer,

we are pleased that you have chosen a product from our range.

Please read these operating instructions carefully and keep them accessible for other users. Please pay attention to the application instructions and notes.

With kind regards

Your Novidion Team

Intended use

- Scope of application: Measurement of blood pressure in humans.
- No contraindications

Main performance and features

- Features:

- Small size
- Light weight
- Portable
- Same effect as mercury sphygmomanometer



- Main Performance Indicators:

- Measuring range: 0 ~ 300 mmHg
- Sub-division: 2 mmHg
- Accuracy: ± 3 mmHg

Adult blood pressure chart (mmHg)		
Recommended Blood Pressure Index		
Classification	Systolic	Diastolic
Optimal	<120	<80
Normal	120-129	80-84
Normal high	130-139	85-89
High mild	140-159	90-99
High moderate	160-179	100-109
Severe	>180	>110




Structure diagram

- Components: Sphygmomanometer is a mechanical elastic element type, which is mainly composed of Palm gauge and cuff parts.
- The patient should be quiet in a relaxed state, preferably without conversation or activity. It is advisable to wear loose clothing.

Usage and measurement methods


- Before using the sphygmomanometer, check whether the palm gauge, the cuff parts, stethoscope are intact; whether the pointer of gauge is located at the "0" position; Whether the response to pressurization is flexible: can be restored to the "0" position by pressing the deflation button at the back of the gauge to ensure the accuracy.
- Take out gauge and cuff from the box, connect and assemble the plastic air pipe on the back of gauge with the single tube of the cuff
- The patient should sit at the table and meditate for 5 minutes before the measurement.
- Adjust the arm height to the height of the heart. If the height is insufficient, place soft pillows under the arm.
- Flatten the arm band so that the ARTERY MARK printed on the band is in the middle.

- Align with the brachial artery, the lower edge of the armband reaches 2 to 3 cm above the elbow socket, and attach it to the upper

arm. It should not be too tight or too loose, and should be able to insert two fingers between the armband and the arm to avoid affecting the correctness.

- The head of stethoscope should be inserted into the cuff and placed on the brachial artery at the elbow.
- Slowly inflate the cuff bladder by squeezing the bulb on the bottom of the gauge. Listen carefully until no brachial artery pulsation can be heard in the stethoscope (The dial of the sphygmomanometer shows a systolic blood pressure about 20-30 mmHg higher than the expected systolic blood pressure), then stop inflation.
- After the pressurization stops, hold the middle of the sphygmomanometer steadily in your hand and gently press the deflation button located at the back of the gauge with the belly of your index finger (the rate of deflation is positively correlated with the strength of pressing the button of the deflating valve), and then slowly release the air in the cuff at a rate of 2 ~ 4 mmHg per second, so that the pointer on the dial will gradually fall at an even speed.
- In the process of descending, the first pulse sound heard through the stethoscope (the thud sound transmitted by the blood flow to reopen the blood vessels and vibrate with them) is recorded as "high pressure" means systolic blood pressure.
- When the pulsatile sound turns into a murmur after a gradual increase in intensity, the tone suddenly becomes muffled and gradually disappears, and when the sound is not audible, it is recorded as "low pressure" means diastolic blood pressure.
- After the diastolic pressure is measured, the deflation valve button can be pressed firmly to accelerate the deflation rate. After the bladder is completely deflated, remove the cuff and complete the measurement process. After the first measurement is completed, deflate completely and wait for at least one minute before repeating the measurement. Generally, the measurement should be repeated 2 to 3 times, and the average of the readings obtained will be the final blood pressure value.

Attention

- Blood pressure in the human body changes throughout the day, so it should be measured at the same time every day.
- After attaching the cuff to the subject's upper arm, the head of stethoscope should be placed in the position of the brachial artery. It is best to hear the pulse sound first, and then inflate and compress it to prevent the pulse sound from not being heard during measurement.
- The measured part of the arm should be at the same level as the heart. If it is lower than the heart, it may result in a blood pressure reading that is higher than the actual reading, otherwise it may be lower than the actual reading.

Maintenance

- In order to maintain accurate readings on the sphygmomanometer, it is necessary to avoid large external impact.
- After using, the sphygmomanometer should be returned to its package and cleaned regularly. It can be wiped with a soft cloth and neutral detergent, and stored in a dry place. Do not fold the rubber hose for too long, fold the cuff into small pieces, or place it next to the stove to prevent excessive temperature from affecting the measure.
- If you have doubts about the accuracy of blood pressure, you can use an accurate mercury sphygmomanometer for calibration. The method is to remove the three-palm catheter on the stethoscope, connect the mercury sphygmomanometer and our gauge at both ends, and connect the cuff at the other end. Using the same pressure to verify whether the readings of the mercury sphygmomanometer and the Palm gauge are the same. If they are different and the values differ by more than 3 mmHg, it means that the Palm gauge is inaccurate. If the readings are the same or the values are within 3 mmHg, only if the pointer has a deviation in position, it will not affect actual use.
- If it exceeds its position, should be calibrated.
- If discover that the pointer cannot back to "0", pls do not adjust the screw without authorization to avoid damaging the movement.
- If a malfunction is found, pls do not open it without authorization. It should be returned to our company or designated service department for repair.
- Storage: Humidity: no more than 80%, Temperature: -5°C + 40°C, no corrosive gases, well-ventilated indoor environment.

Symbol	Meaning	Symbol	Meaning
	Caution! Please refer to the accompanying document (user manual)		Caution! Please refer to the accompanying document (user manual)
	Recyclable		Batch number
	Expiry date		Store in a dry place
	Storage temperature limit		Restriction on storage humidity
	Manufacturer		European representative
	Date of manufacture		Serial number
	Medical device		This product complies with the Medical Device Directives (EU) 93/42/EEC and 2007/47/EEC.

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Novidion GmbH
Lotus NL B.V.
Wuxi Exanovo Medical Instrument Co., Ltd.