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Dear Customer,

we are pleased that you have chosen a product from our range. Please read these operating instructions carefully and carefully and keep them accessible for other users. Please pay attention to the application instructions and notes.

With kind regards
 Your Novidion Team

Precautions

- Do not attempt to repair the Oximeter unless you are professional engineers. Only professionals with maintenance qualification are allowed to perform interior maintenance as necessary.
- Change the contact position between the Oximeter probe and the finger periodically if you are monitoring your SpO₂ levels and pulse rate for a long time (no more than 2 hours).
- Stop immediately if you have broken skin or the blood circulation of your finger is affected during prolong use.
- This product is not designed to be used by newborn babies. Seek for medical care if the measured value goes beyond the normal range and you are sure that the instrument is not malfunctioning.
- The pulse oximeter uses infrared light (invisible to your eyes) to measure your SpO₂ levels. Hence, please do not stare at the light-emitting components of the Oximeter, as that could cause harm and/or potentially blind your eyes.
- This pulse oximeter is not intended to diagnose or treat any medical condition or disease. People who need SpO₂ and pulse rate measurements because of a medical condition should not use the oximeter and should consult with their physician.
- For details about clinical limitations and contraindications, please carefully consult relevant medical literatures.

The following factors may affect the accuracy of the measurement:

- The Oximeter is used in an environment involving high-frequency devices, such as high-frequency electric knives and CT apparatuses.
- Ambient light intensity is too bright. Hence, please avoid direct exposure to strong light (such as beams from operating lamps or sunlight) during measurement.
- The probe of the Oximeter is placed on the same arm that a blood pressure cuff arterial duct or intravenous injection.

- The user suffers from hypotension, severe vascular atrophy, severe anemia, or low oxygen.
- The user is in sudden cardiac arrest or shock state.
- The user is wearing nail polish or artificial nails.
- Do not combine old and new batteries, different brands batteries for using

Warnings

Warning: Do not use the Oximeter in an environment with any flammable gases, flammable anesthetic, or other flammable substances.

Warning: Keep unit and lanyard away from children as the included lanyard may present an entanglement or choking hazard to small children. Adult supervision required; never leave children unattended with unit or lanyard

Warning: Do not throw the batteries into fire, as that could cause an explosion.

Warning: Do not attempt to charge the included batteries, as that could cause leakage, fire disaster, or even explosion. Dispose the used batteries in accordance to the local laws and regulations.


Warning: Do not use the Oximeter in an MRI or CT environment.

Warning: Caution: Do not operate the Oximeter if it is wet. Avoid moving the oximeter from a cold to a hot and humid environment.

Warning: Install the batteries properly before powering on the Oximeter for normal use. Please remove the batteries if you are not planning to use the Oximeter for a long time.

Warning: Close the battery cover when the instrument is in use.

Symbols

Symbol	Meaning	Symbol	Meaning
	Applied part type BF		Date of manufacture
	Caution! Please read this manual.		Manufacturer
% SpO ₂	Oxygen saturation		Authorised representative in the European Community
bpmPR	Pulse rate		Serial number
	No SpO ₂ alarms		Authorised representative in UK
	Follow the instructions for use.		Authorised representative in Switzerland
IP22	The degree of protection against the ingress of water and fine dust		When the end user abandons this product, he/she must send it to the collection point for recycling.

Overview

Oxygen saturation is the percentage of oxyhemoglobin (HbO₂) that is combined with oxygen against all combinable hemoglobin (Hb). It is an important physiological parameter involved in respiration and circulation. The oxygen saturation of arterial blood in a normal human body is 98%. Oxygen saturation is an important indicator of the oxygen condition in the human body. In general, the normal values of oxygen saturation shall not be lower than 94%. If the measured value of oxygen saturation is lower than 94%, an insufficient supply of oxygen is considered.

The pulse rate is the number of pulse beats per minute. Normally, the pulse rate is consistent with the heart rate. In general, the pulse rate of every people is 60 to 90 beats per minute.


The Perfusion Index (PI) usually reflects the limb perfusion status of an examined patient, and shows the detection precision of the instrument as well; that is, examination can still be performed even in the low or weak perfusion condition. The PI of a normal human body is 3% or greater.

Working Principles and Usage

Based on full digital technology, the Finger Pulse Oximeter non-invasively measures the actual content (oxygen saturation) of oxyhemoglobin (HbO₂) in arterial blood using the optical transmittance method.


The Finger Pulse Oximeter measures the blood oxygen saturation and pulse rate of a human body via finger artery. It is applicable to a wide range of fields, such as families, clinics. Use this instrument for measurement before or after sports. You are not advised to use this instrument during sports activities. Do not use it for continuous care for patients.

Schematic Diagram of Display

Schematic Diagram of Display

The following figure shows the information display on the LED screen of the Oximeter in normal detection state:

Power-On button/Function Button Operations

Press the power-on/function button to turn on the oximeter. Once the it is turned on. Simply press or hold the button to perform corresponding

operations.

Press: Press the button for less than 0.5 seconds.

Hold: Press the button for more than 0.5 seconds.

Brightness Setting

Hold the power-on button while the oximeter is in powered-on state, then the oximeter shows a brightness setting interface (as "Interface 1" below shows, "br" represents brightness). Hold the button to adjust brightness. There are 3 brightness settings (1,2,3). 3 is the brightest.

Interface 1

Alert Setting


After setting the brightness, press the power-on button to enter the alert setting interface (as "Interface 2" below shows, "AL" represents alert). Then hold the button to set "AL" to on or off. When "AL" is set to on and the measured values of the blood oxygen saturation and pulse rate go beyond the upper limit or lower limit, the oximeter will beep to alert.


Interface 2

Alert Range Setting


When "AL" is set to on, you can set the upper limit and lower limit of SpO₂ Alert and PR Alert. Press it to switch an option (SpO₂ upper limit, SpO₂ lower limit, PR upper limit and PR lower limit). Hold the power-on button to adjust the limits. As "Interface 3,4,5,6" below show, "Hi" represents upper limit, "Lo" represents lower limit.


Interface 3

Interface 4


Interface 5

Interface 6




Operation Guide

Stick one finger completely into the finger chamber of the oximeter. The fingernail should be facing upward. Release the clip and press the power-on button to power on the pulse oximeter.

⚠ If you do not insert your finger completely into the chamber, measurement will be inaccurate.

⚠ To keep your finger still during measurement. It is also not advisable to use this instrument during sports activities as movement may lead to inaccuracies. Once the reading stabilizes, read the measured values of oxygen saturation and pulse rate on the screen.

NOTE: The oximeter will automatically shut down 10 seconds after you remove your finger.

- ▲ Replace the batteries when the batteries run out of power and the symbol  flickers on the screen.

Install the two AAA dry batteries into the battery slot according to polarity indication, and mount the battery cover.

Cleaning

Power off the instrument and remove the batteries before cleaning. Ensure that the appearance of the instrument is neat, dust-free, and dirt-free. Clean the outer surface of the instrument (including the LED screen) using a piece of dry soft cloth dipped with 75% medical alcohol

Caution: Avoid liquid flowing into the instrument during cleaning.

Caution: Do not immerse any part of the instrument into any liquid.

Disinfection

Before measurement with the instrument, wipe the rubber finger pad using a piece of dry soft cloth dipped with 75% medical alcohol. Clean the finger to be measured using the medical alcohol for disinfection purposes before and after use.

- ▲ Do not disinfect the instrument by means of high-temperature/high-pressure or gas disinfection.

Maintenance

- Remove the batteries from the battery slot and properly store them if you do not plan to use the Oximeter for a long period of time.
- Avoid using the Oximeter in an environment with inflammable gases or using it in an environment where the temperature or humidity is excessively high or low.
- Check the accuracy of the oxygen saturation and pulse rate readings by using an appropriate calibration apparatus.

Technical Specifications

1. Dimensions: 52.6 mm (Length) × 34.1 mm (Width) × 31.2 mm (Height)
Weight: 48.4 g (including two AAA dry batteries)
2. Peak wavelength range of the light emitted from the probe: red light 663 nm ± 3; infrared light 900 nm ± 7.
3. Maximum optical output power of the probe: 60 mW for infrared light (905 nm).
4. Working power supply and maximum current: DC 3.0V; 45mA
5. Normal working condition

Working Temperature	5°C to 40°C (41°F to 104°F)
Relative Humidity	15% to 80%, non-condensing
Atmospheric Pressure	70 kPa to 106 kPa
Rated Voltage	DC 3.0 V

6. Default values and conditions of alert

Parameter	Value
-----------	-------

Oxygen saturation	Upper limit: 100 Lower limit: 94
Pulse rate	Upper limit: 130 Lower limit: 50
Alert condition	When the alert switch is on and the actual measured value goes beyond the preset alert parameter range, the Oximeter gives an alert sound.

7. Technical parameters

Parameter		Value
Display range	Oxygen saturation	35% to 99%
	Pulse rate	25 bpm to 250 bpm
Resolution	Oxygen saturation	1%
	Pulse rate	1 bpm
Measurement precision	Oxygen saturation	±2% (70% to 100%) No requirement (≤ 69%)
	Pulse rate	±2 bpm
Alert range	Oxygen saturation	Upper limit: 50% to 100% Lower limit: 50% to 100%
	Pulse rate	Upper limit: 25 bpm to 250 bpm Lower limit: 25 bpm to 250 bpm
Alert error	Oxygen saturation	± 1% of the preset value
	Pulse rate	The greater of ±10% of the preset value and ±5 bpm

Safety Type

Anti-electric-shock type: internal power supply device

Anti-electric-shock degree: Type BF applied part

Running mode: continuous working

Waterproof grade: IP22

Storage and Transportation

Temperature: -10°C - 50°C (14°F-122°F)

Relative humidity: 10%-93% (no condensation)

Atmospheric pressure: 50kPa-106 kPa.

Product Accessories

- One strap
- Two AAA batteries
- Multilingual operating instructions

Novidion GmbH
 Fuggerstr. 30 • 51149 Cologne • Allemagne
 Tel.: +49 (0) 2203 / 9885 200
 Fax: +49 (0) 2203 / 9885 206
 www.pulox.de • Mail: info@pulox.de

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EMC Information-Guidance and Manufacture's Declaration

1* WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally."

2* WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation."

3* WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ME equipment, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result."

Table 1

declaration - electromagnetic emission	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Not applicable
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable

Table 2

declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Not applicable
Surge IEC 61000-4-5	± 0.5kV, ± 1 kV line(s) to lines ± 0.5kV, ± 1 kV, ± 2 kV line(s) to earth	Not applicable
Voltage dips, short interruptions and voltage variations	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°and 315° 0 % UT; 1 cycle and	Not applicable

on power supply input lines IEC 61000-4-11	70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycles	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m
NOTE: UT is the a.c. mains voltage prior to application of the test level.		

Table 3

declaration - electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Conducted RF IEC 61000-4-6	3 V 0.15 MHz to 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz	Not applicable
Radiated RF IEC 61000-4-3	10V/m 80 MHz to 2.7 GHz	10V/m

Table 4

declaration - IMMUNITY to proximity fields from RF wireless communications equipment					
Immunity test	IEC60601 test level				Compliance level
	Test frequency	Modulation	Maximum power	Immunity level	
Radiated RF IEC 61000-4-3	385 MHz	**Pulse Modulation: 18Hz	1.8W	27 V/m	27 V/m
	450 MHz	*FM+ 5Hz deviation: 1kHz sine	2 W	28 V/m	28 V/m
	710 MHz 745 MHz 780 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m
	810 MHz 870 MHz 930 MHz	**Pulse Modulation: 18Hz	2 W	28 V/m	28 V/m
	1720 MHz 1845 MHz 1970 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m
	2450 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	28 V/m
	5240 MHz 5500 MHz 5785 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	9 V/m

Note* - As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

Note** - The carrier shall be modulated using a 50 % duty cycle square wave signal.

Importeur/ Importer / Importatore / Importador / Importör:

Novidion GmbH, Fuggerstr. 30, 51149 Köln, Germany

Tel.: 02203 – 9885 200, Fax: 02203 – 9885 206

info@pulox.de, www.pulox.de

Shenzhen Jumper Medical Equipment Co., Ltd

ADD: D Building, No. 71, Xintian Road, Fuyong Street, Baoan, Shenzhen, Guangdong, China

E-Mail: info@jumper-medical.com


Tel: +86-755-26692192, 26696279

Web: www.jumper-medical.com

MedPath GmbH

ADD: Mies-van-der-Rohe-Strasse 8,
80807 Munich, Germany

Swiss AR Services GmbH

ADD: Industriestrasse 47, 6300 Zug, Switzerland

E-mail: info@swissarservices.ch

SUNGO Certification Company Limited

ADD: 70 Gracechurch Street,
EC3V 0HR London, United Kingdom

WEEE-Registrierungsnummer: DE24355330

Sollten Sie beim Lesen dieser Anleitung Probleme haben, so können Sie diese auch auf unserer Website downloaden / If you have problems reading this manual, you can also download it from our website / Si vous avez des difficultés à lire ce guide, vous pouvez également le télécharger sur notre site web / Se avete problemi a leggere queste istruzioni, potete anche scaricarle dal nostro sito web / Si tiene

problemas para leer estas instrucciones, también puede descargarlas de nuestro sitio web / Als u problemen hebt met het lezen van deze instructies, kunt u ze ook downloaden van onze website / Jeśli masz problemy z odczytaniem tej instrukcji, możesz ją również pobrać z naszej strony internetowej / Om du har problem med att läsa dessa instruktioner kan du också ladda ner dem från vår webbplats: **<http://www.pulox.de/Anleitungen-Datenblaetter-Retouren>**

Model: JPD-500H(LED)
Version: 1.02 // 15.02.2023