D Digitales Fieberthermometer SC 17

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GB Clinical thermometer SC 17

Special features

The SC 17 is a maximum thermometer, i.e. the highest temperature is displayed after only 60 seconds. The most recent measured value is always saved.

Safety information

- The SC 17 is suitable only for measuring human body temperatures.
- The thermometer contains small parts (batteries, etc) which could be swallowed by small children. Please keep the thermometer and its batteries out of the reach of small children. If the battery or other parts are swallowed you should seek medical help immediately.
- Protect the SC 17 from direct sunshine and temperatures above 50°C.
- Never drop the SC 17, do not bend or open it (with the exception of the battery compartment for battery replacement).

Operating manual

Press the I/O button to switch on. A short signal tone shows that the device is in the "switched on" condition. A segment test takes place simultaneously (checks all display segments and carries out an internal function test (see under "Function test"). After this, if the temperature measured by the temperature sensor is less than 32°C, the display shows an "L" and "°C". After around one second, the "°C" starts to blink. The "L" stands for low (less than 32°C). The current temperature is continuously displayed during measurement, and the "°C" blinks during the entire measurement time. The measurement is completed after around 60 seconds. If the measurement is taken orally or rectally, the approximate final temperature is shown after about 60 seconds and changes only marginaly ater ms time.

After the measurement is completed, the temperature value measured remains in the display until the SC 17 is switched off. To turn the device off after use, press the I/O button briefly. The signal tone will again be sounded. If the device is not switched off, a signal tone is sounded after around 3 minutes and the device then switches off automatically.

Display of last measured value

If you keep the I/O button pressed for longer than three seconds, the last measured value will be shown in the display. As soon as you release the I/O button, the device is ready to take measurements.

Function test

The thermometer is checked for correct functioning automatically every time it is switched on. All segments are displayed, and then "L.C". If a failure (measurement inaccuracy) is recognised, "ERR" (Error) is shown in the display.

Measured value deviation

The measured temperatures will always deviate, depending on the measurement location. This is normal, and caused by deviations in temperature distribution on the human body. You could, for example, measure a rectal temperature of 36,6°C, an oral temperature of 36,1°C and an underarm value of 35,6°C.

Measurement locations

In the mouth (oral)

Small temperature variations in the mouth between measurements are normal. These are caused by varying temperature distribution within the mouth. Try to place the thermometer in one of the two heat pockets under the tongue, to the left and right of the tongue root. The measurement sensor must have good contact with the tissue. The subject should close the mouth and breathe normally through the nose so that the measured result is not affected by breathed air.

Measuring period: approx. 1 minute

In the anus (rectal)

This is the safest method. It is particularly suitable for infants and small children. Guide the thermometer

put carefully max. 2 cm into the anus. Measuring period: approx. 1 minute.

Under the arm (axillary)

From a medical viewpoint this type of measurement leads to inaccurate values. It is therefore unsuitable for carrying out exact measurements.

Measuring the basal body temperature

The device can also be comfortably and practically used for the measurement of fertile or unfertile days in the female cycle using the temperature measurement method.

Changing the battery

As soon as the 🔋 symbol appears in the display, the battery has run out and requires replacement.

Always replace the discharged battery with one of the same type. To do this, proceed as follows:

Pull the cover off the battery compartment. The battery can now be accessed. Remove the use battery with the help of a pointed object (pincers, small screwdriver etc.) and put the new battery (AG 3, LR 41 or SR41) in. Please make sure it is the right way round. You can purchase batteries from jewellers or electrical shops.

Measurement checks

Checking the accuracy of the thermometer is not necessary in most cases. The Medical Product User regulations stipulate a measurement check every two years for devices used in medical practices. If you require a measurement check (at your expense), please send the thermometer to importers address. Please make sure it is well-packaged and that the correct postage stamps are applied.

Cleaning and disinfection

Please only use an alcohol pad. Alcohol pads can be purchased from any dispensing chemist. We recommend the use of Mucocit[®]-F2000 or Helpur[®]-Hplus to disinfect the thermometer when it is being used in medical facilities. Both are effective against all hospital-typical pathogens and can be purchased from any dispensing chemist.

Microbial contamination

Because a hundred percent disinfection is not possible under conditions at home, it is urgently recommended that every patient should use their own fever thermometer.

Guarantee

We give a guarantee period of 2 years from the date of purchase. We guarantee perfect design of this article. We will rectify all material or manufacturing errors during the guarantee period. If, contrary to expectations, defects are found, please send the device back to address of the importer. Please enclose the proof of purchase and a description of the defect so that the problem can be rectified as soon as possible.

Disposal note

Defect devices or devices which are taken out of operation at the end of their service life may still contain residues which could present a risk to the environment. Please dispose of the device in a suitable and correct way at a collection point for old electrical items. Batteries must be disposed of correctly. There are containers for battery disposal at battery dealers and in communal collection points.

Specifications	
Type	maximum thermometer
Measuring range	32,0°C to 43,9°C, Temperature below 32,0°C: display "L" for low. Temperature above 43,9°C: display "H" for high. Automatic internal check on test value 37,0°C.
Measurement accuracy	0,1°C between 35,5°C - 42,0°C at a ambient temperature of 18°C-28°C.
Self-test	Automatic internal check of test value 37,0°C.
Display	In case of deviation > 0,1°C, "Err" will be displayed.
Battery type	3-digit liquid-crystal (LCD display).
Battery type	1 no.: LR 41 or SR 41.
Battery display	The 🔋 shows when the battery is discharged.
Storage temperature	-10°C to 50°C

Explanation of symbols

🏠 Manufacturer	K-jump Health Co., Ltd. No. 56, Wu Kung 5th Rd. New Taipei Industrial Park New Taipei City, 24890 Taiwan
CE0197 CE-sign of a notified body (TÜV Rheinland)	ECREP Polygreen Germany GmbH Ruhlsdorfer Str. 95 14532 Stahnsdorf/Germany
ECREP EU representative	Importeur SCALA Electronic GmbH Ruhlsdorfer Str. 95 14532 Stahnsdorf
🚶 Product type BF	
📅 Disposal notes	
📅 Product period (month/year)	
📖 Read the instruction manual	

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