BE-5000

Lungentrainer - Bedienungsanleitung
Lung trainer - Instruction manual
Entrenador pulmonar - Manual de instrucciones
Allenatore polmonare - Manuale di istruzioni
Entraîneur pulmonaire - Mode d'emploi
Longtrainer - Gebruiksaanwijzing
Lungtränare - Bruksanvisning

www.pulox.de

DE
Novidion GmbH
Fuggerstr. 30 • 51149 Köln • Deutschland
Tel.: +49 (0) 2203 9885-200
Fax: +49 (0) 2203 9885-206
Web: www.pulox.de • Mail: info@pulox.de

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Operating Instructions (EN)

Dear customer,

we are pleased that you have decided to purchase a product from our range. Please read these operating instructions carefully and keep them accessible for other users. Please pay attention to the instructions.

With kind regards

Your Novidion Team

Instructions for use for the Pulox BE-5000 lung trainer

Advantages of lung training

Breathing or lung training serves to strengthen the respiratory muscles and improve lung function with regard to inspiratory breathing (inhalation), which can increase lung volume.

The lung trainer is used, among other things, to help patients who have undergone abdominal or thoracic surgery to improve the depth and duration of their breathing, improve the respiratory function of their recovering lungs, and reduce and prevent post-operative complications.

Its use can also help to increase performance through improved oxygen supply and reduce shortness of breath.

Working principle

When inhaling, resistance is generated inside the lung trainer, which makes breathing in more difficult and thus trains the respiratory muscles.

Features and material:

The lung trainer consists of a housing, an indicator ball, a movable piston, a telescopic tube and a mouthpiece. The D-shaped housing is made of polystyrene, while the telescopic tube, mouthpiece, indicator ball and movable piston are made of polyethylene.

Application:

1. Open the packaging and take out the product.
 2. Remove the telescopic tube with mouthpiece from the handle, pull the telescopic tube to the required length and connect one end of the telescopic tube to the inlet of the lung trainer under the smiley display.
 3. Exhale deeply and let out as much air as possible, then hold the mouthpiece firmly, enclose it with your lips and inhale deeply and evenly so that the indicator ball in the left chamber remains as close as possible to the middle smiley position. The movable piston in the right chamber should rise slowly and remain at the top for as long as possible.
 4. Remove the mouthpiece from your mouth and exhale (don't exhale into device). Repeat steps 3 and 4 for 10 to 15 minutes, then pause and breathe normally.
- Users can use the spirometer for deep breathing training and record the maximum achievable result; patients can perform deep breathing training after surgery under the guidance of medical staff.

Contraindications: No side effects have been identified.

Application notes:

- (1) The numerical value displayed on the lung trainer indicates the suction speed required to display the rise of the ball. For example, '5000 ml' means that the suction speed required to raise the indicator ball is 5000 ml per second. When the indicator ball reaches the top, the maximum suction speed is 5000 ml. The product of the maximum speed and the duration of the indicator ball's ascent correspond to the deep suction volume.
- (2) Use the breathing trainer to perform deep breathing exercises and record the maximum achievable result.

(3) First set the target value for each day, then train for a certain period of time at a low flow rate (e.g. more than 2 seconds, which may take several days depending on lung function), and then, after a certain period of time, increase the inhalation speed for inhalation training until a higher level is reached.

(4) Care and maintenance: Clean the mouthpiece of the product with clean water after each use, dry it well and put it back in the box.

(5) Do not use the product if it has expired or the packaging is damaged.

(6) This product is mainly used for deep training and breathing tests and is not a professional spirometer. The scale on the product surface is a reference value and is for guidance only.

(7) If you feel dizzy or uncomfortable, reduce the intensity of the exercises. If you experience severe discomfort, stop using the product immediately and discuss further use with your doctor.

Storage: Store in rooms with non-corrosive gases and good ventilation.

Symbol	Meaning	Symbol	Meaning
	Follow the instructions for use		Batch number
	Authorized representative in the European Community		Medical device
	This product complies with Directive (EU) 2017/745 on medical devices, a directive of the European Economic Community.		Manufacturer

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**Ningbo Yingmed Medical Instruments Co., Ltd.**
**Lotus NL B.V.**