## GEBRAUCHSINFORMATION

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## INSTRUCTIONS FOR USE

**Aktimed TAPE PLUS** is an elastic self-adhesive tape bandage (kinesiology tape) with Arnica montana 10<sup>6</sup>, Rhus toxicodendron 10<sup>-8</sup> and Ruta graveolens 10<sup>-6</sup> integrated into the adhesive mass. The focus is on the physical and mechanical properties of adhesive tape.

### APPLICATION AREAS

- Kinesiological Taping

### APPLICATION INSTRUCTIONS

- The product can only be used once and needs to be applied by a trained person only.
- Prior to application the skin should be cleaned and fat-free.
- Remove the protection film.
- Apply tape with and without strain.
- The glue unfolds the adhesive effect through body temperature. For this purpose it is advised to smooth out the tape in longitudinal direction from base to end (approx. 60 sec.).
- Apply the tape 30 min prior to sports activity.
- After showering, taking a bath or swimming the tape should be dried with a hair dryer (lowest setting) or by dabbing.
- The tape can occasionally cause skin irritation. In that case please remove the tape.
- The tape can best be removed after taking a shower or a bath.

### SAFETY INSTRUCTIONS

- Do not use in case of an allergy against one of the ingredients, or test the tape on a small area of skin first!
- Do not apply on open wounds, injured skin or unhealed scars!
- Do not use in case of acute or chronic skin diseases (e.g. neurodermatitis, thrush, psoriasis, fungal infection)!
- Do not use in case of Thrombosis or Thrombophlebitis!
- Do not use in case of a peripheral arterial occlusive disease distal to the occlusion!
- Do not use in case of malignant tumors, radiotherapy or chemotherapy!
- If pregnant, please apply the tape only according to your doctor's advice!
- Keep out of reach of children!

### DURATION OF USE

The tape lasts up to 6 days on the skin.

### TECHNICAL INFO

#### Size 5 cm x 5 m:

- AT 20101 PZN 12407103 Aktimed TAPE PLUS light blue
- AT 20102 PZN 12407126 Aktimed TAPE PLUS pink
- AT 20103 PZN 12407155 Aktimed TAPE PLUS black
- AT 20104 PZN 12407161 Aktimed TAPE PLUS beige
- AT 20105 PZN 12407184 Aktimed TAPE PLUS green
- AT 20106 PZN 12407178 Aktimed TAPE PLUS dark blue

**Supporting material:** 95% cotton, 5% spandex

**Adhesive composition:** synthetic rubber based adhesive mass with Arnica montana 10<sup>6</sup>, Rhus toxicodendron 10<sup>-8</sup>, Ruta graveolens 10<sup>-6</sup>

Store in a cool, dry place away from direct sunlight between 8°C and 25°C.

**CE:** Compliant EU Guideline 93/42/EWG for medical products  
**Manufacturer:** Aktimed GmbH, Theaterstr. 18,  
69117 Heidelberg, Germany  
www.aktimed-tape.de

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